

Covidien Valleylab FX8 Electrosurgical Unit – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Covidien Valleylab FX8 Electrosurgical Unit represents the latest advancement in the Valleylab energy family, featuring intelligent TissueFect technology that reads tissue composition in real time. The system monitors and adjusts energy delivery based on tissue impedance for precise performance while ensuring consistent clinical tissue effects using lower energy levels.

The FX8 supports a wide range of surgical procedures requiring division, hemostasis, and dissection. It offers cut, coag, and soft coag modes for surgical flexibility, with Valleylab mode combining monopolar hemostasis and division in a single setting. Shared Coag mode allows two surgeons to operate simultaneously using one energy platform.

The system delivers precise monopolar and bipolar energy with five times better power control than the ForceTriad system and tissue monitoring 2,600 times faster than the Force FX generator, all while maintaining appropriate thermal spread.

FEATURES

- Full Features of an ESU with the LigaSure Vessel Sealing System
- 5 Modes: Cut, Valleylab, Coag, Bipolar, and LigaSure
- 3 LCD Screens, all Touchscreens
- TissueFect™ Sensing Technology
- REM Ground Pad Contact Monitoring
- 4 Coag Modes
- Monopolar and Ligasure Footswitches Included
- Completely Refurbished

SPECIFICATIONS

- Cut: 660 Hz
- ValleyLab Mode: 800 Hz
- Coag: 940 Hz
- Bipolar: 940 Hz
- LigaSure: 440 Hz
- 18"W x 20"D x 10"H
- 30 Lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)