

# Covidien RapidVac Smoke Evacuator – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

The Covidien Valleylab RapidVac Smoke Evacuator efficiently removes odor, particles, and contaminants during electrosurgical procedures. It works with most electrosurgical units, making it a versatile addition to any surgical setup. Surgeons value its quiet operation and ease of use, improving overall surgical workflow and safety in the operating room.

The included SEA3700 25Hr Filter delivers up to 25 hours of high-performance smoke filtration. The RapidVac features a 3-port design that fits various tubing sizes (1/4", 3/8", and 7/8") ensuring flexibility. The 4-stage ULPA filter captures 99.9995% of particles as small as 0.12 microns, creating a clean surgical environment.

Key features include three fan speeds, a filter life meter, and 30% more airflow than previous models. The system includes a power cord, user's guide, pneumatic footswitch (FSSEP), and new ULPA filter (SEA3700), offering powerful, cost-effective smoke evacuation.

## FEATURES

- High Efficiency: Dual-headed pumps ensure a deep vacuum is quickly reached and maintained.
- Quiet Operation: Equipped with a super quiet muffler and insulation system.
- Durability: Features oil-less piston vacuum motors for long-lasting performance.
- Versatility: Suitable for various medical applications beyond liposuction.
- Safety: Designed to prevent back pressure and ensure patient safety.
- Ease of Use: Includes a foot switch for hands-free operation.
- Reliability: If one pump fails, the other three continue to function independently.

## SPECIFICATIONS

- Voltage: 110 or 220 V.
- Weight: 64 lb.
- Dimensions: 17" x 18" x 18".
- Max Vacuum: 29+ inches Hg.
- Max Flow: 7.5 CFM, which is approximately 212.38 LPM.
- Sound Level: 29 ± 3 db.
- Pump: 2 x 2 Piston.
- Horsepower: 2/3 hp.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)