

GE Carescape B450 Patient Monitor with Co2 - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The GE Carescape B450 Patient Monitor is a lightweight, portable monitoring solution weighing only 11 pounds, designed for use across multiple care settings including OR, PACU, NICU, Emergency Care, and Critical Care. Its adaptable design allows easy configuration for different environments, improving workflow and reducing setup time while supporting efficient intra-hospital transport.

The B450 provides essential monitoring parameters including ECG, SpO2, NIBP, CO2, IBP, and temperature measurement. Clinicians can tailor the system to specific patient populations and integrate with EMR systems through optional wireless connectivity and HL7 protocol support, improving data accuracy and care coordination.

Built to handle changing patient acuity levels, the monitor features configurable software to standardize monitoring across departments. Its transport-friendly design ensures continuous patient monitoring during intra-hospital transfers, providing clinicians with reliable tools for safe, effective care delivery in both high-acuity and emergency settings.

FEATURES

- Functions: NIBP, Temp, Pulse Rate and SpO2
- Built-In Pole Mount
- Long run-time rechargeable battery, 8 to 11 hours
- Rugged to withstand heavy use and cleaning
- Multi-Patient Monitor Capabilities
- 6 Months Parts and Labor Warranty
- Adult Cuff (New)
- New Battery

SPECIFICATIONS

- Functions: NIBP, Temp, Pulse Rate, SpO2
- Pole Mount: Built-In
- Capabilities: Multi-Patient
- Battery: Rechargeable
- Battery Run-Time: 8 to 11 hours
- Memory Storage
- 6 Month Parts and Labor Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)