

Hausted HSS500 Surgical Stool – New



Stock Photo - Options May Vary

FACTORY NEW

The Hausted HSS500 Surgical Stool is built for eye, ENT, and plastic procedures. Foot pedals control the hydraulics, so you can adjust the height from 21 to 28 inches without breaking sterility.

In addition, the three-spindle base rides on four dual Tente casters for greater stability, while a front locking system holds the stool in place. For long procedures, the dual-density foam seat with a Dura-Gel insert keeps the user comfortable.

Meanwhile, the Hausted HSS500 Stool's urethane armrest rotates a full 360 degrees, and the armrest support rotates 180 degrees for stable arm positioning. Additionally, the contoured backrest adjusts to fit each user, and the ABS thermoformed hood makes cleaning easy.

FEATURES

- **Gas Lift Adjustment:** Allows for quick, one-touch height changes to suit different clinicians
- **360-Degree Swivel:** Provides a full range of motion for easy movement around patients
- **High-Impact Base:** Features a five-leg design for maximum stability during procedures
- **Durable Vinyl Upholstery:** Resists stains and stands up to frequent medical-grade disinfectants
- **Compact Footprint:** Fits easily into tight exam rooms or small laboratory stations

SPECIFICATIONS

- **Weight Capacity:** 250 lbs
- **Seat Height Range:** 16.5" - 21.75" (standard gas lift)
- **Seat Diameter:** 16"
- **Base Diameter:** 22"
- **Base Material:** High-impact black plastic

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)