

# Hausted Horizon Series – Replacement Stretcher Pad – 30" x 77" x 4"



*Stock Photo – Options May Vary*

## FACTORY NEW

Hausted Horizon Series – Replacement Stretcher Pad – 30" x 77" x 4"- The standard stretcher pads are made by combining three layers of high density, high resilient, medical grade urethane foam helps provide outstanding pressure redistribution, postural support and comfort. Take patient care to the next level with these durable stretcher pads. If you have any questions live chat with us or give us a call we are always happy to help.

This pad for the Hausted Horizon Series has the following measurements: 30" Width 77" Length 4" Thickness This pad is available in different options. Call today to place your order.

Note! – Loop Velcro down center bottom unless otherwise specified on order.

## FEATURES

- Universal Corners: 2 Round / 2 Angled (Foot-Radius Cut, Head-angle cut)
- Fluid Protection Zipper Guard ensures protection against fluid contamination.
- Anti-microbial, Anti-static, Stain resistant, Abrasion resistant, and Non-allergenic.
- Fire rated CA 117
- Other options available

## SPECIFICATIONS

- Width: 30"
- Length: 77"
- Thickness: 4"
- Weight: 25 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)