

Skytron 6702 Standard Series Replacement Pad Set



**For use with the Skytron 6702
Stock Photo - Options May Vary*

FACTORY NEW

Skytron 6702 Standard Series Replacement Pad Set. The Skytron 6702 Standard Series Replacement Set suits low risk patients undergoing short duration procedures. This pad set also fits Skytron models 6700B, 6701, and 6702.

This cushion set uses medium firm, high density foam to deliver steady, reliable support throughout routine surgical cases. The conductive vinyl cover resists fluids and static while meeting O.R.

safety standards. Facilities choose this standard set when they do not need advanced pressure relief.

FEATURES

- Latex-Free
- Designed for Short-Duration Procedures
- Velcro Kit Complete with Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant

SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)