

Skytron 6702 Premium Series Replacement Pad Set



**For use with the Skytron 6702
Stock Photo - Options May Vary*

FACTORY NEW

Skytron 6702 Premium Series Replacement Pad Set. The Skytron 6701 Premium Series Replacement Set delivers lasting patient comfort during longer duration procedures. This pad set also fits Skytron models 6700B, 6701, and the 6702.

This cushion set combines high density foam with Visco memory foam, reducing pressure points and supporting proper patient positioning. The conductive Fusion III vinyl cover stretches four ways and resists fluids, tears, and hospital disinfectants.

Surgical teams rely on this premium set to maintain comfort and performance during extended cases.

FEATURES

- Designed for Longer-Duration Procedures
- Velcro Kit Comes Complete with the Pad
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable

SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)