

Skytron 6302 Standard Series Replacement Pad Set



**For use with the Skytron 6300/6302
Stock Photo - Options May Vary*

FACTORY NEW

The Skytron 6302 Standard Series Replacement Pad Set suits low risk patients undergoing short duration procedures. This cushion set uses medium firm, high density foam to deliver steady, reliable support throughout routine surgical cases.

The conductive vinyl cover resists fluids and static while meeting O.R. safety standards.

Facilities choose this standard set when they do not need advanced pressure relief. Standard Series Pad Set Features: 2" thick overall Constructed with medium-firm high-density foam Covered with electrically conductive vinyl Color: Black If you need padding for Bariatric patients or longer cases, or just added patient comfort, take a look at our premium pad for this table here.

FEATURES

- Designed for Short-Duration Procedures
- Hook & Loop Kit comes included with the Pad Set
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable

SPECIFICATIONS

- Construction: Medium Firm High-Density Foam
- Exterior Cover: Electrically Conductive 2-Way Stretch 'Lectrolite Vinyl
- Pad Thickness: 2"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)