

# Skytron 6302 Premium Series Replacement Pad Set



\*For use with the Skytron 6300/6302  
*Stock Photo - Options May Vary*

## FACTORY NEW

The Skytron 6302 Premium Series Replacement Pad Set delivers lasting patient comfort during longer duration procedures. This cushion set combines high density foam with Visco memory foam, reducing pressure points and supporting proper patient positioning.

Furthermore, the conductive Fusion III vinyl cover stretches four ways and resists fluids, tears, and hospital disinfectants. Surgical teams rely on this premium set to maintain comfort and performance during extended cases.

## FEATURES

- For Longer Duration Procedures
- Waterfall Edge Construction
- Fluid and Stain Resistant
- X-ray Permeable
- Anti-Microbial

## SPECIFICATIONS

- Construction: Medium Firm High-Density with 1" Visco memory foam
- Cover: Conductive 4 Way Stretch Fusion III Soft Vinyl
- Pad Thickness: 2"

\*Brochure and images are a general overview of the product. Actual product options may vary.

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)