

Verathon BVI 9400 Bladder Scanner – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Verathon BVI 9400 Bladder Scanner – Refurbished. The BVI 9400 is a portable 3D ultrasound that quickly measures bladder volume and post-void residual (PVR), without surgery. It works on both adults and small children, so one device covers a wide range of patients.

Additionally, NeuralHarmonics technology makes the scanner accurate and easy to use, and it works without a sonographer. Using the BVI 9400 is simple. With one click, the device captures ultrasound images from multiple angles inside the body.

Then it builds those images into a 3D dataset of the bladder, and calculates the volume right on screen.

FEATURES

- Portable, 3D ultrasound bladder volume measurement
- NeuralHarmonics technology for accurate results
- No sonographer required
- Works on adults and small children
- One-click image capture and volume calculation

SPECIFICATIONS

- Dimensions: 48 inch cord length
- Display Type: LCD Display
- Power Source: Battery Operated
- Probe Frequency: 1-3 MHz
- UNSPSC Code: 42201711

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)