

# Allen Stirrup Cart – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

Allen Stirrup Cart – Refurbished. This cart gives you a strong, compact way to store and protect your stirrups between procedures.

It holds stirrups or leg holders from any manufacturer, not just Allen, so it fits whatever your OR already has. The molded, four caster polymer base is dent resistant, so it holds up well to daily bumps and rolling between rooms.

Because the cart is compact, it tucks easily into a storage room or hallway without taking up much space. A deep storage bin keeps clamps and pads organized and secure, so nothing gets lost or damaged between uses.

## FEATURES

- Strong, compact, lightweight cart
- Holds and transports stirrups or leg holders from any manufacturer
- Dent resistant polymer base
- Deep storage bin for clamps and pads
- Two hook sizes for a range of stirrup rod diameters

## SPECIFICATIONS

- OEM Part Number: 931101
- Base: Molded four caster polymer, dent resistant
- Small Hook: Fits rods up to 3/4" (2 cm)
- Large Hook: Fits rods up to 1" (2.5 cm)
- Lead Time: 5 to 7 business days

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)