

# Kendall SCD 700 Sequential Compression System – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

Kendall SCD 700 Sequential Compression System, Refurbished. This controller helps prevent deep vein thrombosis (DVT) and pulmonary embolism (PE) with sequential, gradient compression to the leg, foot, or both at once.

It uses Vascular Refill Detection technology to read each patient's own vascular refill time. The system then adjusts therapy to match, instead of running one fixed cycle for everyone.

This refurbished unit includes a new battery and a new tubing set, so it's ready to put into service right away. Please note, sleeves come as a separate purchase.

## FEATURES

- Sequential, gradient, circumferential compression
- Vascular Refill Detection (VRD) technology
- Patient detection and animated alerts
- IP23 rated, Xylex housing, protected battery
- Quiet operation with vibration dampeners

## SPECIFICATIONS

- Manufacturer: Covidien (trade name: Kendall)
- Condition: Refurbished
- Includes: New Battery, New Tubing Set
- Sleeves: Sold separately
- Ingress Rating: IP23

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)