

Verathon GlideScope AVL Video Laryngoscope – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Verathon GlideScope AVL Video Laryngoscope – Refurbished. This system is a trusted, reusable video laryngoscope built for high-volume ORs, ICUs, and trauma centers.

It pairs the reusable GlideScope GVL blade, angled at 60 degrees, with a high resolution color monitor. Furthermore, because the monitor uses Reveal anti-fog technology and a built-in LED light, you get a clear, wide-angle view that helps cut down on failed intubations, even in difficult airways.

The AVL system runs on a digital platform, unlike the standard GVL system's analog platform, so you also get real-time video recording and snapshot capture built in.

FEATURES

- 60 degree angulated reusable GVL blade
- Reveal anti-fog technology with integrated LED
- Real-time video recording and snapshot capture
- Export to external displays
- Rugged, reusable design for daily use

SPECIFICATIONS

- Manufacturer: Verathon
- Condition: Refurbished
- Dimensions: 6.8 in H x 8.8 in W x 3.1 in D (17.4 x 22.3 x 8 cm)
- Weight: 2.2 lbs (1 kg)
- Display: TFT Color VGA, 640 x 480 px, 6.4 in (16.3 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)