

Welch Allyn Connex 6000 Series Vital Signs Monitor – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Welch Allyn Connex 6000 Series Vital Signs Monitor – Refurbished. This vital signs monitor, also known as a patient monitor, delivers accurate bedside readings straight to your EMR. The refurbished unit features SureBP 15-second NIBP measurement, SureTemp Plus or Braun ThermoScan thermometry, and Masimo or Nellcor SpO2.

Additionally, a large color touchscreen simplifies training and daily operation on busy care floors. Built for adult, pediatric, and neonatal patients, this vital signs device supports interval monitoring with customizable alarm programs. Furthermore, clinicians can manually enter weight, height, respiration rate, and pain scores.

Customizable Early Warning Scoring, including MEWS and PEWS, enables clinical decision support right at the bedside.

FEATURES

- Comfortable for your patients
- Practical for the team: One device for all vital signs & standardized supplies
- Ease of use to optimize your time with the patient
- Intuitive menu structure
- Single-button operation

SPECIFICATIONS

- Completely Refurbished
- Patient Range: Neonatal, Pediatric, & Adult
- Functions: NIBP, Temp, Pulse Rate, and SpO2
- Built-in Pole Mount
- 6 Months Parts and Labor Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)