

Mayfield Horseshoe Headrest Pad – New



Stock Photo – Options May Vary

FACTORY NEW

Mayfield Horseshoe Headrest Pad: NH9002. The Mayfield Adult Horseshoe Headrest Pad Set, also known as a horseshoe headrest cushion or cranial positioning pad, is a new, non-OEM direct replacement for the Mayfield Horseshoe Headrest. This complete set includes both left and right pads, restoring cushioning and patient contact points without the wait or cost of OEM parts.

It's built to fit seamlessly onto your existing headrest frame. Compatible model of headrest include the A-1011 and A-1012. Horseshoe headrests support patients in the prone or supine position during neurosurgery, spine surgery, and other delicate procedures.

Fresh padding helps maintain proper weight distribution and reduces pressure points during long cases.

FEATURES

- Complete pad set, includes both left and right pads
- New, non-OEM direct replacement for the Mayfield Horseshoe Headrest
- Designed for prone and supine patient positioning
- Supports neurosurgery, spine surgery, and other cranial stabilization procedures
- Fits existing Mayfield Horseshoe Headrest frames without modification

SPECIFICATIONS

- Product Number: NH9002
- Condition: New
- Configuration: Complete set (left and right)
- Compatibility: Mayfield Horseshoe Headrest A1011, A1012

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)