

Stryker Concave Eye Stretcher Headrest Pad Replacement – New



Stock Photo – Options May Vary

FACTORY NEW

Stryker Concave Eye Stretcher Headrest Pad Replacement. The Stryker Concave Eye Stretcher Replacement Pad is a new, direct-fit pad for Stryker 1068, 1069, 1079, and 1089 eye surgery stretchers. The concave shape cradles the patient's head securely during ophthalmic procedures, keeping it centered and stable throughout delicate surgical positioning.

Velcro attachment makes installation and removal simple, while the contoured design helps reduce pressure points during longer eye surgery cases. Eye surgery stretchers require specialized headrest support that standard flat pads can't provide.

This pad's concave contour accommodates the patient's head naturally, improving comfort and stability during ophthalmic procedures.

FEATURES

- New concave headrest pad for Stryker eye surgery stretchers
- Fits Stryker models 1068, 1069, 1079, and 1089
- Velcro attachment for secure, stable positioning
- Available in 2", 3", and 4" thicknesses
- Black color only

SPECIFICATIONS

- Compatibility: Stryker 1068, 1069, 1079, 1089
- Color: Black
- Attachment: Velcro
- Condition: New
- Lead Time: 7-10 business days

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)