

Precision Medical 3300 PM Series 3-Mode Continuous-Intermittent Vacuum Regulator



Stock Photo - Options May Vary

FACTORY NEW

Precision Medical 3300 PM Series 3-Mode Continuous-Intermittent Vacuum Regulator. The Precision Medical 3300 PM Series regulator provides reliable, accurate cycling on and off during intermittent suctioning. Additionally, this series has a peak vacuum limit feature to prevent adjusting the vacuum too high.

Furthermore, it's ideal for sensitive applications where limiting max vacuum is critically important. Moreover, the intermittent timing cycle integrated into the Precision Medical 3300 is factory set at approximately 16 seconds on and 8 seconds off. As a result, staff gain consistent, predictable suction cycling for demanding procedures.

Specifically, fully adjustable, 2-to-1 ratio timing cycles allow further customization for specific clinical needs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)