

Venaflow Elite Aircast SCD System – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Venaflow Elite Aircast SCD System – Refurbished. The VenaFlow Elite System is designed as a prophylaxis for deep vein thrombosis. Additionally, it combines two technologies, rapid inflation and graduated sequential compression, to mimic ambulation and accelerate venous velocity.

Furthermore, breathable calf, foot, and thigh cuffs enhance patient comfort. Also, these cuffs work with one universal pump for easy setup. Moreover, this system meets ACCP and AAOS guidelines for hospital and surgery center use.

As a result, facilities gain a trusted, evidence-backed DVT prevention option. Specifically, automatic cuff detection identifies and adjusts pressure based on the attached cuff configuration. In addition, preset pressures and alarms eliminate the need for manual adjustments.

FEATURES

- Automatic Cuff Detection with Preset Pressures
- Telescoping Bed Hanger (Extends up to 3.5")
- Universal Pump for Calf, Thigh, and Foot Cuffs
- Compliance Alarm for Patient Monitoring
- Low-Profile, Lightweight Design

SPECIFICATIONS

- Indication: Prophylaxis for Deep Vein Thrombosis (DVT)
- Models: 30B (Standard), 30B-S (Slow Fill), 30B-SB (Slow Fill + Battery Backup)
- Bed Hanger Extension: Up to 3.5"
- Compliance: Meets ACCP and AAOS Guidelines
- Battery-Installed Units: Available Upon Request

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)