

Champion 8500 Treatment Table

FACTORY NEW

Champion 8500 Treatment Table. The Champion 8500 Treatment Table features a durable dove-gray powder-coated steel frame built for daily clinical use. Additionally, this table includes a paper holder with cutter, adjustable leveling feet, and a comfortable patient surface.



Stock Photo - Options May Vary

Furthermore, standard height is 30 inches, with specialty heights of 19, 25, or 36 inches also available. Moreover, the 2-inch medium-density foam surface is upholstered in health-grade vinyl for easy cleaning. As a result, staff can quickly wipe down this table between patients.

Specifically, the Champion 8500 table supports a 400 lb weight capacity, suitable for a wide range of patients. In addition, adjustable leveling feet help keep the table stable on uneven floors.

FEATURES

- Features an adjustable backrest, two storage compartments, two sliding doors, and an adjustable shelf in each compartment
- Easy-clean laminate surfaces
- Ships assembled
- Seamless rounded corner top
- Heavyweight, knit-backed upholstery

SPECIFICATIONS

- 9044-24 Dimensions: 72" L x 24" D x 31" H
- 9044-27 Dimensions: 72" L x 27" D x 31" H
- 9044-30 Dimensions: 72" L x 30" D x 31" H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)