

Champion 8570 Adjustable Treatment Table

FACTORY NEW

Champion 8570 Adjustable Treatment Table. The Champion 8570 Adjustable Treatment Table includes a powder-coated steel frame, paper holder, and adjustable leveling feet. Additionally, this table features comfortable 2-inch medium-density foam upholstered with health-grade vinyl.



Stock Photo - Options May Vary

Furthermore, a gas spring-operated headrest uses a one-touch control for quick, easy adjustment. Also, this design helps staff position patients faster during exams. Moreover, standard height is 30 inches, with specialty heights of 25 or 36 inches also available.

As a result, the Champion 8570 table fits a wide range of clinical spaces and staff needs. Specifically, the headrest angle adjusts from 6 to 57 degrees for added flexibility. In addition, adjustable leveling feet keep the table steady on almost any floor.

FEATURES

- Features an adjustable backrest, two storage compartments, two sliding doors, and an adjustable shelf in each compartment
- Easy-clean laminate surfaces
- Ships assembled
- Seamless rounded corner top
- Heavyweight, knit-backed upholstery

SPECIFICATIONS

- 9044-24 Dimensions: 72" L x 24" D x 31" H
- 9044-27 Dimensions: 72" L x 27" D x 31" H
- 9044-30 Dimensions: 72" L x 30" D x 31" H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)