

Tuttnauer T-Classic 9-Liter Portable Autoclave



Stock Photo - Options May Vary

FACTORY NEW

Tuttnauer T-Classic 9-Liter Portable Autoclave – New. This compact, portable steam autoclave fits well in small clinics or as a backup for larger practices. It brings Tuttnauer performance with a simple, easy-to-use design.

In addition, it runs fully automatic cycles with one-touch operation, starting at just 18 minutes. At only 11.7 pounds, the T-Classic 9 is the lightest tabletop autoclave on the market. Meanwhile, it needs very little counter space, which makes it easy to fit almost anywhere.

An electronic detector shuts off power if water runs low. A built-in thermal fuse also guards against overheating and electrical faults.

FEATURES

- A simple one-step operated device
- With One Touch, the T-Classic 9 will begin an automatic cycle until completion
- An electronic detector will turn off the power if the water level is insufficient
- Perfect autoclave for limited space
- The T-Classic 9 offers competitively quick unwrapped cycle times, in as short as 18 minutes*

SPECIFICATIONS

- Max Width / Height (overall): 13.4" / 13"
- Sterilizing temperature: 250°F
- Max Net Weight: 5.3 kg/11.7 lbs.
- Min sterilizing time: 18 min.
- Chamber Diameter/Height: 8.2"/9" Voltage: 120V

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)