

Tuttnauer EZ9 Fully Automatic Autoclave – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

Tuttnauer EZ9 Fully Automatic Autoclave – Refurbished. This autoclave makes sterilization easy with one-touch operation. It fills, sterilizes, exhausts, and dries on its own with a single button.

In addition, it works well for clinics, dental offices, and medical practices that need reliable, hands-free sterilization. The Refurbished Tuttnauer EZ9 Fully Automatic Autoclave features a closed-door active drying system that keeps items sterile while drying both packs and pouches fully. Meanwhile, a HEPA 0.2 μ m replaceable air filter boosts safety and cleanliness by cleaning air entering the chamber during drying.

Staff can also set and save cycle settings to match their exact needs.

FEATURES

- Closed-door active drying system to maintain sterility and ensure efficient drying of packs and pouches
- A High-efficiency air pump circulates hot air after the sterilization cycle for fast drying
- Hepa filter 0.2 μ m air filter (replaceable)
- Pre-loaded cycles for unwrapped instruments, wrapped/packs, glassware, and additional drying cycle
- All program parameters can be customized and stored

SPECIFICATIONS

- Chamber Dimensions: 9" x 18"
- Chamber Volume: 5 Gallons
- Overall dimensions: 21.5"D x 20"W x 14.4"D
- Standard Cassettes Capacity: 2 Full and 2 Half
- Tray Dimensions: 16.3"D x 6.7"W x 0.8"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)