

# Tuttnauer T-Top 11 Autoclave



Stock Photo – Options May Vary

## FACTORY NEW

Tuttnauer T-Top 11 Autoclave – New. This autoclave serves small to medium dental and medical clinics. It covers the basic needs of general clinic sterilization while helping boost daily productivity.

In addition, this fully automatic tabletop unit offers strong capacity in an affordable, reliable package. The Class S model uses gravity steam sterilization, while Class B uses a pre/post vacuum method that keeps air from being exhausted during the cycle. Meanwhile, Class S units can upgrade to Class B through a simple USB software update, so clinics don't need to buy a new device.

This unit is plug-and-play, needing only a power source to run.

## FEATURES

- Fast Cycles: The Fast Cycle times for wrapped instruments, complete cycle, and including drying.
- Optimal Chamber Capacity: More than 75% chamber capacity with 5 extra-long trays.
- Convenient loading: Convenient and easy loading with door opening to 110°.
- Advanced documentation: A variety of documentation and traceability options and barcodes.
- User-friendly touchpad control panel: Displays programs, cycle times, water levels, and more.

## SPECIFICATIONS

- Volume: 6 Gallon/23 L
- Overall Dimensions (DxWxH): T-Edge 10: 21.1" x 19" x 19.7"
- Chamber Size (Depth x Ø): T-Edge 10: 18" x 10" diameter
- Standard Cassette Capacity: T-Edge 10: 3 Full + 3 Half
- Tray Capacity: 5 Trays

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)