

Protekt 100 Series Replacement Foam Mattress



Stock Photo - Options May Vary

FACTORY NEW

Protekt 100 Single Layer Pressure Redistribution Foam Mattress - New. This mattress gives superior pressure redistribution for low- to moderate-risk patients. Its foam core uses a high-resilience 1.5 density formulation at a 28 ILD rating, allowing deep immersion and even weight support.

In addition, this mattress carries Medicare code E0184. The nylon cover is vapor permeable, which cuts down on perspiration and gives a low-shear, low-friction surface. Meanwhile, this cover is also anti-microbial, anti-bacterial, and fluid resistant.

A 360-degree zipper makes the cover easy to remove, clean, or replace, and an anti-contamination flap keeps fluids away from the foam core.

FEATURES

- Superior pressure redistribution for low to moderate risk patients
- Deep immersion foam core for even weight distribution
- Low shear, low friction cover surface
- Easy-clean, removable cover with anti-contamination flap
- Optional raised rails available on all sizes

SPECIFICATIONS

- Foam: High-resilient 1.5 density, 28 ILD
- Cover: Vapor-permeable, anti-microbial, anti-bacterial nylon, fluid resistant
- Closure: 360-degree zipper with anti-contamination flap
- Fire Rating: Meets Federal Fire Code 16 CFR Part 1633
- Weight Capacity: 350 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)