


Protekt 300 Series Replacement Foam Mattress



Stock Photo - Options May Vary

FACTORY NEW

Protekt 300 Replacement Foam Mattress. This mattress gives excellent pressure redistribution for moderate- to high-risk patients. Its multi-tiered design uses a top layer of high-density, high-resilience die-cut foam in the torso and head sections.

In addition, this mattress carries Medicare code E0184. The sloped heel section pairs with open-cell memory foam to reduce pressure and boost blood flow by shifting load from the heel to the lower leg. Meanwhile, the PTSS rail system gives extra support for easy patient transfer and edge-of-bed sitting.

Furthermore, a stretch nylon top cover reduces perspiration with a low-shear, low-friction surface.

FEATURES

- Superior pressure redistribution for low to moderate risk patients
- Deep immersion foam core for even weight distribution
- Low shear, low friction cover surface
- Easy-clean, removable cover with anti-contamination flap
- Optional raised rails available on all sizes

SPECIFICATIONS

- Foam: High-resilient 1.5 density, 28 ILD
- Cover: Vapor-permeable, anti-microbial, anti-bacterial nylon, fluid resistant
- Closure: 360-degree zipper with anti-contamination flap
- Fire Rating: Meets Federal Fire Code 16 CFR Part 1633
- Weight Capacity: 350 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)