

Protekt 500 Series Replacement Hospital Bed Mattress



Stock Photo - Options May Vary

FACTORY NEW

Protekt 500 Pressure Redistribution Hospital Bed Foam Mattress - New. This mattress uses VGPT (Visco Gel Polymer Technology) for pressure redistribution well beyond traditional foam. It also moves heat away from the body much faster, cutting down on perspiration.

In addition, its open-cell design gives a real boost in breathability compared to standard foam mattresses. This design of the Protekt 500 offers more support and less cradling at higher compression, spreading pressure over a wider surface area. Meanwhile, a 30-degree heel slope protects the vulnerable heel area by shifting pressure to the leg.

The PTSS rail system also gives extra support for easy patient transfer and edge-of-bed sitting.

FEATURES

- Rapid heat dissipation, reduces perspiration and pressure sore risk
- Superior pressure redistribution at higher compression levels
- 30-degree heel slope to protect vulnerable heel areas
- PTSS rail system for easy transfer and edge-of-bed sitting
- Optional raised rails available on all sizes

SPECIFICATIONS

- MPN: 81051 (36"x76"x6"), 81052 (36"x80"x6"), 81053 (36"x84"x6")
- Foam: VGPT (Visco Gel Polymer Technology) open-cell gel-infused foam
- Cover: Stretch nylon top, anti-microbial, anti-bacterial, fluid resistant
- Bottom Cover: Triple laminate vinyl, fluid resistant, fire retardant
- Fire Rating: Meets Federal Fire Code 16 CFR Part 1633

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)