

Oakworks Shield | 82611

FACTORY NEW

The Oakworks Shield is designed to reduce the risk of contamination via aerosols and droplets by providing a barrier between caregivers and patients. It is also effective in separating patients from one another in larger rooms.

Portable and easy to move, the Oakworks Shield can be rolled into place when needed. Adjustable in height, you can easily find the balance between access and safety. The optically clear film can be easily positioned to create a barrier between two people while not restricting access to the patient.

The shield provides months of crystal clear protection from a single roll. Instead of expensive permanent installations that become cloudy over time, the shield is flexible to meet the needs of caregivers performing exams and procedures. Shield film and support elements can be disinfected and reused, minimizing cost.



Stock Photo - Options May Vary

FEATURES

- The flexible shield creates an effective barrier to aerosols and droplets, even when proximity is a problem
- The shield can be used between patients to minimize exposure to aerosols and droplets. The shield is also easily moved when not needed
- To create a new barrier, just cut off the old barrier, pull out the new film, and clip it in place
- Shield film and support elements can be disinfected and used over and over again, minimizing cost
- The shield's optically clear film is very flexible, tear, and abrasion-resistant
- The unit comes partially disassembled-assembles in a few minutes

SPECIFICATIONS

- Barrier Size 48" x 36" Shield Area
- Height Range 8" - 38" (bottom of the shield)
- Base 5 Star with Wheels
- Shield Support Arm Carbon Fiber with 2 Clips
- Support Pole Stainless Steel 2 Piece Pole
- Replacement Film Rolls Roll makes 25 barriers, Part # 82664
- Weight 11 lbs. without Roll, 22 lbs. with Roll
- Warranty 1 Year

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)