

Harloff SureDry Tambour Door 14 Scope Drying Cabinet w/Key Lock | SC8036TKDP-14



Stock Photo - Options May Vary

FACTORY NEW

The Harloff SureDry Tambour Door 14 Scope Drying Cabinet delivers powerful, efficient drying for medical scopes. It uses a positive-pressure HEPA-filtered drying system that ensures clean, dry air circulation. Two high-volume intake fans and one exhaust fan continuously replace internal air, promoting rapid, hygienic drying.

This reliable system helps protect valuable scopes from moisture damage and contamination. Designed for modern healthcare environments, the SureDry cabinet combines performance with convenience. The roll-top tambour door saves space and allows quick access to all stored scopes.

The key lock adds secure storage, keeping equipment safe between uses. Hospitals and clinics trust Harloff for durable, easy-to-maintain solutions.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)