

Amico Five-Function Electrical Bed | AC-400C



Stock Photo - Options May Vary

FACTORY NEW

Amico Five-Function Electrical Bed: AC-400C. Amico Corporation, founded in 1974, designs, manufactures and markets a broad range of products for every health care facility department. Furthermore, their mission is to provide state-of-the-art products that facilitate the caregiver's role in providing the best outcomes for patients.

The Amico Apollo Hospital Bed Collection gives caregivers the control they need in the recovery environment. The simple design and distinctive features bring efficiency to today's medical and surgical settings. The Amico Apollo AC 400C features easily removable head and footboards that give caregivers easy access in an emergency.

The head and footboards feature an exclusive polypropylene blend for uncomplicated cleanup and durability.

FEATURES

- Easily removable head and footboards
- Instant CPR release lever
- 550 lb (249.48 kg) weight capacity
- LED night light
- Bed bumpers for added protection

SPECIFICATIONS

- Model Number: AC- 400C
- Overall Length: 88" (224 cm)
- Overall Width: 42" (107 cm)
- Patient Surface: 35" x 77" (89 x 196 cm)
- Height Range: (to the surface of the litter), High 29" (74 cm), Low 18" (46 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)