

Stryker InTouch ICU Hospital Bed – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Stryker InTouch Critical Care Bed combines intuitive controls with advanced technology to enhance patient safety, mobility, and caregiver efficiency. Intuitive controls on all sides, a powerful touch screen interface, and one-touch electric brake put full control at your fingertips. The bed supports early mobility with a 16" low bed height, independent knee gatch, and one-touch cardiac chair positioning.

Siderails encourage proper transfer techniques and safe positioning. Features include protocol reminders, sound therapy, and spoken language translations to ensure comfort and communication throughout care.

The bed features iBed Wireless connectivity for seamless integration with hospital information systems, enabling fall prevention, asset tracking, and bedside documentation. The integrated Zoom Drive System improves transport efficiency, while iBed Awareness alerts staff to safety risks. Open architecture allows the bed to adapt to facility needs.

FEATURES

- Touch screen interface
- Smart scale
- Chaperone Bed Exit System
- Electric brake with mechanical back-up
- CPR foot pedal and CPR-ready headboard
- Enhanced cardiac chair
- One-handed siderail release
- Three-zone bed exit alarm

SPECIFICATIONS

- 90" (228.6 cm)
- Siderails up 42" (106.7 cm)
- Siderails down 40.5" (102.9 cm)
- 550 lb. (249 kg)
- Low (standard with 6" casters) 16" (40.6 cm)
- Backrest 0°–70°
- Knee gatch 0°–15°
- 35" x 84" (89 x 213 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)