

# Stryker S3 3005 MedSurg Hospital Bed - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

Stryker S3 MedSurg Hospital Bed, Model 3005. The S3® bed uses modern technology and smart design to improve patient care and make work easier for caregivers.

This fully electric bed works well in many hospital settings. It features one-handed side rails and simple controls, making position changes easy for staff and comfortable for patients.

Stryker Secure 3 (Model 3005) Features: Key features include a low bed height for safer transfers, a retractable bed frame for comfort during repositioning, and removable head and footboards with integrated control panels. Additionally, these design elements make caregiving more efficient and patient experiences more comfortable.

## FEATURES

- Centrally-located brake pedal
- 4-wheel, steel brake with swivel lock system
- 6" (15 cm) casters
- One-touch 30-degree head of bed angle and lock
- One-hand, intermediate siderails

## SPECIFICATIONS

- Siderails Up: 41.5" (105.4 cm).
- Siderails Down: 39.5" (100.3 cm).
- High: 30" (76.2 cm).
- Low (standard with 6" casters): 16" (40.6 cm).
- Backrest: 0 to 60°

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)