

# UMF Power Exam Chair | 4011P



Stock Photo - Options May Vary

## FACTORY NEW

The UMF 4011P Power Exam Chair is an ultra procedure chair designed for healthcare providers in specialties including OB/GYN and Oncology. It features the OneTouch Patient Positioning® System with both hand and foot control, offering four power functions: hi-lo, back, tilt, and leg lift with extension. The chair meets ADA and U.S. Access Board regulations to accommodate most patients safely.

The chair accommodates a 650 lb. patient weight capacity and provides low access at 18.5" minimum height. It includes a double-articulating headrest, removable seamless upholstery with Prefixx® Protective Finish, stainless steel drain pan, and stirrups with four lateral positions. The design supports multiple clinical positions including Trendelenburg, Lithotomy, and Semi-Fowler.

Manufactured in the U.S.A., the 4011P is engineered to protect patients from infection, facilities from damage, and staff from injury. It comes with a three-year limited warranty and features a hospital-grade electrical receptacle with child safety cover.

## FEATURES

- 650 lb. (294.84kg) patient weight capacity
- Low access ADA height 18.5" ( ±0.5") 46.9cm (±1.27cm)
- four power functions: hi-lo, back, tilt, leg lift & extension
- Hand control with OneTouch Patient Positioning® System, two programmable positions & return-to-home Hand Control
- Foot control with OneTouch Patient Positioning® System, two programmable positions & return-to-home
- Removable, seamless upholstery with Prefixx® Protective Finish
- Double-articulating headrest with adjustable height
- Stainless steel drain pan

## SPECIFICATIONS

- Height: min. 18.5" (± .5") 46.9cm (±1.27cm) at seat max. 37" at the seat (93.9cm)
- Seat width: 28" (71.12cm)
- Back width: 26.5" (67.31cm)
- Length: 73"-80" with foot section extended (185.42cm - 203.2cm)
- Leg rest: 28" W x 11.5" D (71.12cm W x 29.21cm D)
- Stirrups: four lateral positions, adjustable length
- Electrical: 120 VAC, 60HZ
- Net weight: 465 lb (211 kg)

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)