

# Blickman Kay Cabinet | 7953SS



Stock Photo - Options May Vary

## FACTORY NEW

Blickman 7953SS - Kay Cabinet. Blickman's stainless steel cabinets feature a seamless welded face that prevents dirt collection and provides long-lasting durability with unmatched strength and corrosion resistance.

Their freestanding design sits on four legs with glides or can move easily on casters. The Kay (shown) and Kennedy cabinets feature two separate sections with unlocked or locked compartments.

The Kay Cabinet is designed for efficient storage and easy access, with ample space for medical instruments, supplies, and equipment. Its customizable compartments offer flexibility for various storage needs, ensuring organization and security.

## FEATURES

- Freestanding
- Separate upper and lower sections.
- 16" deep upper section:
- One locked glass swinging door on two hinges
- Three glass shelves

## SPECIFICATIONS

- Dimensions: 30" x 58" x 16"
- Weight: 219 lbs
- Construction: 304 stainless steel
- 1/4" clear tempered sheet glass
- Welding: All welding by Heliarc method and exposed welds are ground smooth

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)