

# Clinton Pediatric Series/Aquarium, Blood Drawing Chair | 66011-AQ

## FACTORY NEW

Clinton Pediatric Series/Aquarium, Blood Drawing Chair: 66011-AQ. Clinton Industries has produced quality medical equipment designed for superior clinical performance for over three decades. Clinton is one of the broadest and best sources for a wide range of quality, practical medical equipment.

Furthermore, our medical equipment and accessories lines are easily accessible, coordinated, and cost-effective. Clinton's 66011-AQ Blood Drawing Chair Features: Custom high pressure, Clinton Kids™, Aquarium flip arms, Easy-clean laminate arm surfaces with black, T-molded edges, and Graphics embedded in the laminate.

Cleaning with standard disinfectants does not affect them. It also features a ClintonClean™ short armrest, Padded Contoured Seat, and Thick Padded Backrest.



Stock Photo - Options May Vary

## FEATURES

- Custom high pressure, Clinton Kids™, Aquarium flip arm, and right armrest
- Easy-clean laminate arm surfaces with black, T-molded edges
- Graphics are embedded in the laminate and are not affected by cleaning with standard disinfectants
- ClintonClean™ short armrest
- Thick padded seat and backrest

## SPECIFICATIONS

- Seat (width, depth, height): 17.5" X 18" D X 26" H
- Overall: 31" X 27" X 42"
- Arm Height Range: 33" X 39"
- Arm Adjustments: Six

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)