

Clinton Pediatric Arctic Circle, Blood Drawing Chair | 66011-AC



Stock Photo - Options May Vary

FACTORY NEW

Clinton Pediatric Arctic Circle, Blood Drawing Chair: 66011-AC. For over thirty years, Clinton Industries has produced quality medical equipment designed for superior clinical performance. Clinton is one of the broadest and best sources for a wide range of quality, practical medical equipment.

Our medical equipment and accessories lines are easy to use, coordinated, and cost-effective. Clinton Pediatric Series/Arctic Circle Blood Drawing Chair features custom high pressure, short Clinton Kids™, Arctic Circle flip arm, and right armrest. In addition, the chair is an Easy-clean Laminate arm surface with black, T-molded edges, and graphics are embedded in the laminate and are not affected by cleaning with standard disinfectants.

FEATURES

- Custom high-pressure, Clinton Kids™, Arctic Circle flip arms
- Easy-clean laminate arm surfaces with black, T-molded edges
- Graphics are embedded in the laminate and are not affected by cleaning with standard disinfectants
- ClintonClean™ short armrest
- Thick padded seat and backrest

SPECIFICATIONS

- Model Number: 66011-AC
- Load capacity under normal use: 400 lbs
- Seat 17.5" W x 18" D x 26" H
- Overall Dimensions: 31" W x 18" D x 26" H
- Arm Height Range: 33" — 39"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)