

# Riester ri-focus™ LED Headlight System | 6090

**FACTORY NEW**

Riester ri-focus™ LED Headlight System: 6090.



*Stock Photo - Options May Vary*

This headlight provides ideal illumination of the examination area with a 40 cm distant, focus-able, white and cold light with 5500 Kelvin and 7100 Lux.

Clearly distinguishable color nuances, stark contrasts, as well as high-performance lighting to enable a faster and easier diagnosis. Long-life lithium batteries or NiMH rechargeable batteries plus plug-in charger - no power cord needed during the examination.

## FEATURES

- 40 cm distant, focus-able, white and cold light.
- 5500 Kelvin and 7100 lux
- The LED lamp employed has a life of approximately 15,000 hours
- Can be supplied with power using standard commercial AAA batteries, long-life lithium batteries or NiMH rechargeable batteries plus plug-in charger
- No annoying power cord need to be used during the examination

## SPECIFICATIONS

- White LED with 6,000 Kelvin, 140 lumen
- LED with up to 100,000 Lux
- Uses 2 Lithium batteries size CR 123, 4 batteries size AAA (alkaline), or 4 rechargeable batteries size AAA (NiMH)
- Weight: 210 grams (without batteries), 242 grams with CR 123 or 259 grams with rechargeable AAA

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)