

# Harloff Bronchoscope and ENT Small Scope Drying Cabinet, 10 Capacity, Key Lock | SC5430DRDP



Stock Photo - Options May Vary

## FACTORY NEW

The Harloff SC5430DRDP Bronchoscope and ENT Small Scope Drying Cabinet with 10 capacity and key lock. holds up to ten scopes with an adjustable mounting system.

Furthermore, it uses a positive-pressure HEPA-filtered fan system with two intake and two exhaust fans to dry scopes. The Harloff SC5430DRDP replaces the air inside four times per minute, and the fan run time is adjustable.

The Harloff SC5430DRDP includes the best-in-class Positive Pressure HEPA Filtered Drying System and features a key lock. You can Trust Harloff for durable, high-performance carts built to support your healthcare facility's requirements.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)