

Tuttnauer T-Classic 9-Liter Portable Autoclave



Stock Photo - Options May Vary

FACTORY NEW

Tuttnauer T-Classic 9-Liter Portable Autoclave. The Tuttnauer T-Classic 9 Autoclave is a lightweight, portable, and cost-effective steam sterilizer designed specifically for small to medium-sized medical and dental practices. With a 9-liter capacity, this compact autoclave offers a simple, fully automatic one-touch operation, making sterilization fast, easy, and affordable.

The T-Classic 9 is perfect for small clinics because of its compact size. It starts sterilization with one simple touch, making work easier and reducing mistakes.

A water level detector keeps the device safe by turning it off if the water is too low.

FEATURES

- A simple one-step operated device
- With One Touch, the T-Classic 9 will begin an automatic cycle until completion
- An electronic detector will turn off the power if the water level is insufficient
- Perfect autoclave for limited space
- The T-Classic 9 offers competitively quick unwrapped cycle times, in as short as 18 minutes*

SPECIFICATIONS

- Max Width / Height (overall): 13.4" / 13"
- Sterilizing temperature: 250°F
- Max Net Weight: 5.3 kg/11.7 lbs.
- Min sterilizing time: 18 min.
- Chamber Diameter/Height: 8.2"/9" Voltage: 120V

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)