

Harloff M-Series Baby Scale Cart | MDS3024SC



Stock Photo - Options May Vary

FACTORY NEW

Harloff M-Series Baby Scale Cart: MDS3024SC. This M-Series cart holds weighing scales for infants and babies on top, with ample room below for storing PPE, blankets, and pediatric supplies in the open equipment area.

Additionally, this cart is perfect for NICUs, Pediatric Clinics, and Children's Hospitals! The mobile baby scale cart features 5" Tente casters for easy transport of the infant scale and supplies. Furthermore, it also includes a 3" high drawer and a 21" open storage area with seven adjustable shelf height positions.

Additionally, Other Cart Features Include Ergonomic Side Handles, Protective Bumpers, and Full Extension Drawer Slides.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)