

Tuttnauer T-Edge 10 Autoclave



Stock Photo - Options May Vary

FACTORY NEW

Tuttnauer T-Edge 10 Autoclave. The Tuttnauer T-Edge 10 is the only sterilizer upgradeable from gravity to a pre/post vacuum system. As a result, this feature lets your clinic start with basic sterilization and easily upgrade to advanced Class B performance.

It saves costs and maximizes long-term value for medical and dental practices. The T-Edge 10 meets the latest ANSI/AAMI ST55 sterilization standards for top efficiency. Not only that, it is built strong and reliable.

Furthermore, its design ensures high efficiency and low operating costs, making it perfect for busy clinics that need trusted infection control.

FEATURES

- **Fast Cycles:** The Fast Cycle times for wrapped instruments, complete cycle, and including drying.
- **Optimal Chamber Capacity:** More than 75% chamber capacity with 5 extra-long trays.
- **Convenient loading:** Convenient and easy loading with door opening to 110°.
- **Advanced documentation:** A variety of documentation and traceability options and barcodes.
- **User-friendly touchpad control panel:** Displays programs, cycle times, water levels, and more.

SPECIFICATIONS

- **Volume:** 6 Gallon/23 L
- **Overall Dimensions (DxWxH):** T-Edge 10: 21.1" x 19" x 19.7"
- **Chamber Size (Depth x Ø):** T-Edge 10: 18" x 10" diameter
- **Standard Cassette Capacity:** T-Edge 10: 3 Full + 3 Half
- **Tray Capacity:** 5 Trays

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)