

# Tuttnauer T-Edge 11 Autoclave



Stock Photo - Options May Vary

## FACTORY NEW

Tuttnauer T-Edge 11 Autoclave. The Tuttnauer T-Edge 11 is a unique and innovative sterilization solution designed to adapt to your practice. As the only gravity and airflow removal autoclave on the market that is upgradable to pre/post vacuum (Class B), the T-Edge 11 offers unmatched versatility, allowing you to scale your sterilization capabilities as your clinic grows.

The T-Edge 11 meets ANSI/AAMI ST55 standards, ensuring top safety, strong performance, and long-term reliability. It maximizes efficiency and minimizes operational costs, making it a smart investment for modern dental and medical facilities.

## FEATURES

- **Fast Cycles:** The Fast Cycle times for wrapped instruments, complete cycle, and including drying.
- **Optimal Chamber Capacity:** More than 75% chamber capacity with 5 extra-long trays.
- **Convenient loading:** Convenient and easy loading with door opening to 110°.
- **Advanced documentation:** A variety of documentation and traceability options and barcodes.
- **User-friendly touchpad control panel:** Displays programs, cycle times, water levels, and more.

## SPECIFICATIONS

- **Volume:** 7 Gallon/27 L
- **Overall Dimensions (DxWxH):** T-Edge 11: 21.1" x 20" x 19.7"
- **Chamber Size (Depth x Ø):** T-Edge 11: 18" x 11" diameter
- **Standard Cassette Capacity:** T-Edge 11: 4 Full + 4 Half
- **Tray Capacity:** 5 Trays

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)