

NUVO V800 Surgical Table



Stock Photo - Options May Vary

FACTORY NEW

NUVO V800 Surgical Table. For nearly 20 years, our sole focus at Nuvo has been to provide the medical field with the highest quality, most cost-effective surgical lights, tables, equipment management systems, video integration, and other products for hospital operating rooms.

All of which can greatly help you make a difference in the care of your patients. NUVO V800 Features: an exceptionally sturdy, low-profile base that helps ensure maximum surgeon comfort to reduce fatigue during long procedures.

Additionally, features unprecedented C-arm access, self-diagnosis with an Instant Error Code display, and a control system that automatically prevents sectional conflicts.

FEATURES

- Unmatched versatility. Exceptional quality. Remarkable value
- • 800 lb. Weight Capacity
- • Removable, Split Leg Section
- • Kidney Bridge
- • Sliding Top

SPECIFICATIONS

- Max. Weight Capacity: 800 lbs (363 kg)
- Height: 26.8 - 44 in (68 - 112 cm)
- Slide: 12.2 in (31cm)
- Lateral Tilt: 20°
- Trendelenburg: 30°

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)