

# Skytron 6702 Surgical Table – Refurbished



Stock Photo – Options May Vary

## FULLY REFURBISHED · CLINICALLY CERTIFIED

Skytron 6702 Hercules Surgical Table – Refurbished. The Skytron 6702 Hercules Surgical Table is the strongest surgical table on the market with an impressive 1,200 lb lift capacity. This makes it ideal for bariatric surgeries and complex procedures. Additionally, its 210° top rotation enhances surgical access without repositioning the patient.

This feature improves patient safety and gives every member of the surgical team better access. With its rotating tabletop, the 6702 Hercules simplifies imaging access at both the head and foot ends. Furthermore, the low-profile base stays clear of imaging equipment like the C-arm.

This design reduces setup time and increases surgical efficiency. As a result, operating rooms can handle more specialties without extra effort.

## FEATURES

- 210° top rotation
- 1,200 lb lift and 1,000 lb articulation weight capacity
- Removable back and leg sections simplify accessory setup
- Power beach chair positioner allows for automatic, one-touch setup of shoulder procedures
- Table height adjustable from 23-5/8" to 41-1/2"

## SPECIFICATIONS

- Length w/o headrest 66 3/8", 1,665mm
- Length with headrest 79", 2,006mm
- Width (back/seat plate) 19 5/8", 500mm
- Width with side rails 21 7/8", 555mm
- Height 23 5/8 to 41 1/2", 600-1,053mm

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)