

Steris 4085 Surgical Table – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED · CLINICALLY CERTIFIED

The Steris 4085 Surgical Table is a multi-purpose operating table designed for strength, reliability, and modern functionality across diverse surgical settings. With a 1,100 lb weight capacity in normal position, it accommodates a wide range of procedures. The table features wired hand control, battery operation, hydro-electric brakes, and a sliding tabletop for C-arm compatibility.

Key movements include height adjustment, backrest and leg section positioning, Trendelenburg, lateral tilt, tabletop slide, adjustable headrest, and powered kidney lift. The table includes Return to Level functionality and North American side rails that fit most surgical accessories.

This refurbished unit comes with a new pad set and has been professionally disassembled, mechanically refurbished, and cosmetically restored to ensure reliable performance in clinical environments.

FEATURES

- C-Arm Compatible Table Top
- Wired Hand Control
- Battery Operation
- Hydraulic Floor Locks
- Return to Level
- North American Side Rails Fit Most Accessories
- 1,100 lbs Weight Capacity in Normal Position
- Completely Refurbished

SPECIFICATIONS

- Weight Capacity: 1100 lbs in Normal Position, 500 lbs in Reverse Orientation
- Table Length (Overall): 81"
- Table Width (Overall): 21 5/8"
- Tabletop Slide Range: 9" (227 mm) to head, 9" to foot
- Patient Surface Width: 20"
- Pad Thickness: 2"
- Power Requirements: 120V, 4.5 Amps

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)