


Steris 5085 Surgical Table – Refurbished



Stock Photo – Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The Steris 5085 Surgical Table is designed for maximum versatility, easily adapting to nearly any surgical specialty. With an 18" sliding tabletop and 1,200 lb. capacity, this table eliminates the need for multiple specialty tables. The large imaging area, low height, and removable leg section enhance access and positioning for various surgical procedures.

Built for performance and convenience, the Steris 5085 features a telescoping stainless-steel column for easier cleaning, along with self-leveling floor locks for added stability on uneven surfaces. The lightweight and impact-resistant design ensures long-lasting use in high-demand environments.

The easy-to-read LCD hand control allows quick, precise adjustments with minimal effort during critical procedures. The table supports up to 1,000 lbs. during full articulation and 600 lbs. with articulation and slide, with a height range from 22.8" to 42.5".

FEATURES

- Sliding top of 18"
- Maximum weight capacity of 1,200 lbs. (patient support, raise/lower on centered column)
- Light weight and impact resistant
- Removable leg section with no protruding pins
- Telescoping stainless-steel column for easy cleaning
- Self-leveling floor locks for stability on uneven floors
- Easy-to-read LCD hand control
- Large imaging area

SPECIFICATIONS

- Table length: 81"
- Table width: 20"
- Table height range: 22.8" to 42.5"
- Tabletop slide range: 9" to head & to foot
- Weight capacity (on centered column): 1,200 patient support (lower/raise), 1,000 lbs. (full table articulation), 600 lbs. (full table articulation, slide)
- Trendelenburg/Reverse: 30°
- Lateral tilt (left/right): 20°
- Head section: +/- 90°

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)