

# Getinge 5624 Warming Cabinet - Refurbished



Stock Photo - Options May Vary

## FULLY REFURBISHED • CLINICALLY CERTIFIED

Getinge Warming Cabinet - Model 5624. The Getinge 5624 warming cabinet is an outstanding choice for facilities that need efficient warming of blankets and IV fluids but have limited floor space. Designed to maximize functionality in compact environments, it offers the convenience of a double-door configuration without occupying excessive space.

The Getinge 5624 Warming Cabinet features clear glass door panes that allow staff to quickly check inventory levels at a glance, reducing unnecessary door openings and preserving consistent internal temperatures. This helps improve energy efficiency and maintain optimal warming conditions for patient comfort and safety.

## FEATURES

- "Look-in" glass doors provide visual inventory control
- Large internal volume for storage of solution bottles and blankets
- Forced air convection heating
- LED display will alert operator of an interruption of

## SPECIFICATIONS

- Exterior Dimensions: 30"W x 60"H x 24.7"D
- Interior Dimensions (Upper): 26"W x 13"H x 22.7"D
- Interior Dimensions (Lower): 26"W x 40"H x 22.7"D
- Capacity: 4.1 Cubic Feet
- Capacity: 13.4 Cubic Feet

Toll-Free

Local

800-627-3215 | 813-740-8640

salesteam@venturemedical.com

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)