

# Getinge 5524 Tabletop Warming Cabinet - Refurbished

## FULLY REFURBISHED · CLINICALLY CERTIFIED

Getinge Castle 5524 Warming Cabinet. The Getinge 5524 Warming Cabinet provides a reliable way to warm saline solutions and blankets in patient care settings.



*Stock Photo - Options May Vary*

It's ideal for areas that require quick access and precise temperature control, like surgery, recovery, OB/GYN, ICU, and ER. This compact tabletop unit combines convenience and performance in a space-saving design that fits easily into busy clinical environments.

The Castle 5524 features a sleek design that includes a tempered glass "look-in" door, allowing staff to visually check inventory without opening the cabinet and disrupting internal temperatures. Additionally, the convection-heated chamber ensures even heat distribution, keeping fluids and linens consistently warm and ready for use.

## FEATURES

- Tabletop Model
- Compact Design
- Digital Display and Controls
- Standard 120V Plug
- Stainless Steel Construction

## SPECIFICATIONS

- Exterior Dimensions: 30"W x 17.5"H x 24.7"D
- Interior Dimensions: 26"W x 13.5"H x 22.7"D
- Capacity: 4.1 Cubic Feet
- Power: 115V, 60 Hz
- Current Draw: 4 Amps

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)