

Tuttnauer T-Edge 10B Autoclave



Stock Photo - Options May Vary

FACTORY NEW

Tuttnauer T-Edge 10B Autoclave. The Tuttnauer T-Edge 10B Autoclave is a powerful, compact Class B sterilizer designed to meet the strict sterilization needs of modern dental and medical practices. Engineered for performance and efficiency, the T-Edge 10B features advanced pre- and post-vacuum sterilization cycles, delivering deep steam penetration to thoroughly sterilize even complex, hollow instruments.

Ideal for practices that value speed and reliability, the T-Edge 10B offers fast 36-minute wrapped instrument cycles, including drying. This significantly reduces turnaround time while maintaining high sterilization standards.

Furthermore, over 80% optimal chamber capacity fits longer trays, boosting cycle productivity and helping busy clinics meet high demand.

FEATURES

- Pre-Post Vac - Class B Autoclave
- Larger chamber capacity and shorter cycle times
- Convenient and easy loading with door opening to 110°
- Advanced documentation and traceability options and barcode
- Optimal chamber capacity (23 Liter/6 Gallons)

SPECIFICATIONS

- Chamber Size: 18" x 10"
- Volume: 6 gal (23L)
- Dimensions (L x W x H): 22.8" x 18.9" x 19.23"
- Standard Cassette Capacity: 3 full & 3 half
- Trays: 5

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)