

# Tuttnauer T-Edge 11B Autoclave



Stock Photo - Options May Vary

## FACTORY NEW

The Tuttnauer T-Edge 11B is a Class B benchtop autoclave designed for dental and medical practices requiring fast, reliable sterilization. It features advanced pre- and post-vacuum steam cycles that fully sterilize complex and hollow instruments, completing a full wrapped cycle with drying in just 36 minutes to minimize downtime.

With a 27-liter chamber capacity and 110° door opening, the T-Edge 11B accommodates extra-long trays and provides ergonomic loading and unloading. The unit supports multiple water filling and drainage methods for flexible installation options.

The autoclave includes Wi-Fi connectivity for remote monitoring via smartphone or tablet, a user-friendly touchscreen interface, and advanced documentation features with barcode capability to meet sterilization record-keeping requirements.

## FEATURES

- Pre-Post Vac - Class B Autoclave
- Convenient and easy loading with door opening to 110°
- Advanced documentation and traceability options and barcode
- Remote Wi-Fi monitoring via smartphone or tablet
- Range of automatic one-button programs

## SPECIFICATIONS

- Chamber Size: 18" Depth x 11" Diameter
- Volume: 7 gal (27L)
- Dimensions (L x W x H): 21.1" x 19.68" x 19.68"
- Standard Unwrapped Cycle Time: 20 Minutes
- Voltage: 110V or 220/230V

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)