

GE MAC 5500 HD EKG Machine - Refurbished



Stock Photo - Options May Vary

FULLY REFURBISHED • CLINICALLY CERTIFIED

The GE MAC 5500 HD EKG Machine (Refurbished) offers both comprehensive 12-lead analysis and built-in interpretation capabilities. Patient data and results are displayed on a 10.4" backlit color AM LCD screen. It features a full alphanumeric keyboard for easy entry of patient demographics.

With one-touch controls for printing, saving, and transmitting data, the MAC 5500 HD simplifies operation and improves workflow efficiency. Another advantage is that this MAC 5500 HD supports full-page printouts and employs age- and gender-specific interpretation criteria for accurate 12SL diagnoses. Additionally, the MAC 5500 HD includes QRS and P wave signal averaging to enhance the precision of its measurements.

*An optional rolling cart is available for the MAC 5500 HD.

FEATURES

- 10" Color Display
- Marquette Hookup Advisor
- Marquette 12SL with Gender-Specific interpretation
- Simultaneous 12/15-lead acquisition, storage, and assessment
- Internal Modem

SPECIFICATIONS

- Display Type: 10.4 in (264 mm) diagonal graphics backlit color AM LCD
- Display Resolution: 640 x 480 pixels with waveform enhancement
- Display Data: Heart rate, patient name, ID, clock, waveforms, lead labels, speed, gain and filter settings, warning messages, prompts, and help messages
- Height: 3.7 in (9.4 cm) with display closed
- Width: 15 in (38.1 cm)

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)